

Dynamic Multi-Modality Fused Imaging, Analysis, Computer Aided Diagnosis System




Figure 1

202090-8216001

Multi-Modality Imaging System with Common Focused 2D Curved Detector




Figure 2

202050-0216001

Overall Multi-Modality Imaging System Block Diagram




Figure 3

X-ray & Focused 2D Curved Detector Arrangement




Figure 4

202090-8216200F

APPROVED	O.G. FIG.
BY	CLASS
CRAFTSMAN	SUBCLASS

Cone Beam Source Collimation & Cone Beam Shaped Filter







Figure 5

202507 321600T

X-ray Cone Beam Focal Spot - Curved Detector Optics

Curved Detector to reduce spatial resolution loss and Best Conversion efficiency of X-ray





FIG. 6a

Focal spot from V-groove Type Anode has similar spot size appearance




FIG. 6b




FIG. 6c

Figure 6

2 Dimensional Focal Spot Dithering for Improved Cone Beam

Spatial Resolution

X-ray Focal Spot Geometric Dithering For Doubling the Spatial Sampling Rate




Figure 7a

2D X-ray Focal Spot Dither Pattern for 3D Cone Beam VCT




Figure 7b

Normal
Sampling

Dither to Improve
Sampling


Figure 7b

2 times the sampling rate

Focused 2D
Curved
Detector


Figure 7c

Spatial Resolution comparison between Single Sampling and 2X Dither Sampling



X-ray Optical
System Response
before Sampling

Figure 7e



Normal Nyquist Sample
Freq. & aliased optical
response

Figure 7

202050-62162001

Focused 2D Curved Detector Module

View Showing Focused 2D
Anti-scatter Collimation with
2D Focused Pixels

Focused Curved Detector
Module







FIG. 8b

Figure 8

+

†



Impulse Response Shaping from Rectangular to Variable gaussian Roll-off Function. Shaping may be Fixed or Controlled



Figure 9

Multi-Modality XGA Detector Module

X-Ray Mode



FIG. 10a

PET Mode




FIG. 10b

NM/SPECT Mode




FIG. 10c

Figure 10

Detector Module Multi-Modality Collimation




Figure 11

XGA Detector Module Signal Processing






Figure 12

System with Optional PET Anti-Scatter Baffle





Figure 13

APPROVED	O.G. FIG.
BY	CLASS SUBC.
CRAFTSMAN	

202050-3216E001

PET - Anti-Scatter Baffle SEPTA



Z-Axis 3D Mode Scatter Reduction Baffles




Figure 14

System With Cone Beam Focused NM/SPECT Collimation




Figure 15

APPROVED	O.G. FIG.	
BY	CLASS	SUBCLASS
DATE		

+

20E050" 22F6E001

NM/SPECT Mode with Collimation Ring





Figure 16

+

202050" 32162001

Cone Beam NM/SPECT LEHR Collimation and Focused 2D Curved

Detector Array



X-ray Gamma Ray Area Detector. [XGA] Detector Which is Focused at Point Where X-ray Focal Spot is.

Figure 17

X-Ray Detector Scatter Rejection with Focused 2D Curved Collimation




Figure 19

Sequencing of X-ray Sources for Adaptive Scatter Correction




Figure 20

Modulation and Demodulation for Scatter Correction with Multiple Sources




Figure 21

System Level Diagram of Modulation and Demodulation For Multiple

Sources for VCT




Figure 22

20E090-82F6001

Step and Shoot VCT Imaging




Figure 23

Spiral VCT with Multiple Heads




FIG. 25a




FIG. 25b

Spiral Path with 3 Heads with
respective Central Rays on
Reconstruction Cylinder




Figure 25




FIG. 25c

Detector Projection to Scan Center




FIG. 25d

Virtual Detector Sampling Region

Cone Beam Slant Source Collimation for Spiral VCT Imaging




Figure 26

Multi-Plane Planning System Imaging




Figure 27

Whole Body Dose Control From Planning System




FIG. 28a

Figure 28

FIG. 28b

202050-32156001

Dynamic Timing Control




Figure 29

Retrospective Gated Imaging System




Figure 30

Prospective Gating Control System with Cardiac EKG





Figure 31

Prospective and Retrospective Gated DAQ and Reconstruction Imaging

Prospective Gating Control





FIG. 32a

Retrospective Gating Control



Multiple cardiac cycles to fill needed Views. Collect all views in (n) cycles of Heart.

FIG. 32c




FIG. 32b

Multi Cycle - Contiguous




FIG. 32e

Figure 32

Gated DAQ and Reconstruction for Retrospective Cine' Dynamic Cardiac Imaging




FIG. 33a




FIG. 33b




FIG. 33c

Figure 33


PET Transmission, Attenuation & Scatter Correction

VCT Attenuation MAP




Transmission Attenuation
Map at 511 KEV Energy Level from VCT Images

PET Attenuation Correction



Correction Map for PET New Corrected PET
Projections for OSEM Recon.

PET Scatter Correction



Scatter Correction from VCT Images and
Count Rates on a Projection View Basis

Figure 34

202090-321600F


APPROVED	O.G. FIG.	
BY	CLASS	SUBC
RAFTSMAN		

NM/SPECT Transmission, Attenuation & Scatter Correction


VCT Attenuation MAP



Transmission Attenuation
Map at NM/SPECT Energy Levels from VCT Images



Correction Map for NM/SPECT New Corrected
SPECT Projections for OSEM Recon.



Scatter Correction from VCT Images and
Count Rates on a Projection View Basis

Figure 35

20250326001

Patient Fused Multi-Modality Imaging and Analysis System




Figure 36

Interventional Image Control System




Figure 37

+

2025090" 02T600F

Multi-Modality Imaging with Independent X-Ray VCT, PET, and NM/SPECT Image Acquisition System




Figure 38

+

APPROVED	O.G. FIG.	
BY	CLASS	SUBC
RAFTSMAN		

202050* 8216E001

Multi-Modality Imaging with Independent X-Ray Single Head VCT, PET, and NM/SPECT Image Acquisition System




Figure 39

APPROVED	O.G. FILE
BY	CLASS
CRAFTSMAN	SUBJECT

+

202050-8262001

Multi-Modality Imaging with Independent X-Ray 4th Generation VCT, PET, and NM/SPECT Image Acquisition System




Figure 40

+

APPROVED	O.G. FIG.
BY	CLASS
DRAFTSMAN	SUBJECT

Multi-Modality Imaging System with Stationary
Focused 2D Curved Detector for VCT, PET and NM/SPECT Imaging




Figure 41

20250-826001

†
Multi-Modality Imaging with Common Gantry and Independent X-Ray VCT,
 PET, and NM/SPECT Image Acquisition System




Figure 42

+

202507 8266001

APPROVED	O.G. FIG.
BY	CLASS
DRAFTSMAN	30200

**Multi-Modality Imaging with Common Gantry and Independent X-Ray
Single Head VCT, PET, and NM/SPECT Image Acquisition System**




Figure 43

202050-8216001

APPROVED	O.G. FILE	
BY	CLASS	SUB
CRAFTSMAN		

**Multi-Modality Imaging with Common Gantry and Independent X-Ray
4th Generation VCT, PET, and NM/SPECT Image Acquisition System**




Figure 44

APPROVED	O.G. FIG.
BY	CLASS SUBC.
DRAFTSMAN	

202090-82166001

Multi-Modality Imaging with Common Gantry and Independent Single X-Ray 4th Generation VCT, PET, and NM/SPECT Image Acquisition System




Figure 45